

Product Datasheet

Prepro-Neuromedin S (70-103) (orb372737)

Description	Peptide to Prepro-Neuromedin S (70-103)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Phe-Leu-Phe-His-Tyr-Ser-Arg-Thr-Gln-Glu-Ala-Thr-His-Pro-Val-Lys-Thr-Gly-Phe-Pro-Pro-Val-His-Pro-Leu-Met-His-Leu-Ala-Ala-Lys-Leu-Ala-Asn
Formula	$C_{180}H_{271}N_{49}O_{44}S$
Purity	> 95%
Source	Synthetic
MW	3857.50
Expiration Date	12 months from date of receipt.