

Product Datasheet

Prodynorphin (226-238) (orb372746)

Description	Peptide to Prodynorphin (226-238)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Gly-Gly-Phe-Leu-Arg-Arg-Gln-Phe-Lys-Val-Val-Thr
Formula	$C_{74}H_{115}N_{21}O_{17}$
Purity	> 95%
Source	Synthetic
MW	1570.87
CAS Number	85006-82-2
Expiration Date	12 months from date of receipt.

