

Product Datasheet

pro-Epsilon-Tx1X/12 (orb372748)

Description	Peptide to pro-Epsilon-Tx1X/12
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Leu-Lys-Arg-Thr-Ile-Arg-Thr-Arg-Leu-Asn-Ile-Arg
Formula	$C_{66}H_{126}N_{26}O_{16}$
Purity	> 95%
Source	Synthetic
MW	1539.90
Expiration Date	12 months from date of receipt.