

Product Datasheet

proFIX18 (orb372749)

Description	Peptide to proFIX18
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Val-Phe-Leu-Asp-His-Glu-Asn-Ala-Asn-Lys-Ile-Leu-Asn-Arg-Pro-Lys-Arg
Formula	$C_{95}H_{157}N_{31}O_{27}$
Purity	> 95%
Source	Synthetic
MW	2165.50
Expiration Date	12 months from date of receipt.