

Product Datasheet

proFIX28 (orb372750)

Description	Peptide to proFIX28
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Val-Phe-Leu-Asp-His-Glu-Asn-Ala-Asn-Lys-Ile- Leu-Asn-Arg-Pro-Lys-Arg-Tyr-Asn-Ser-Gly-Lys-Leu- Glu-Glu-Phe-Val
Formula	$C_{149}H_{235}N_{43}O_{44}$
Purity	> 95%
Source	Synthetic
MW	3332.80
Expiration Date	12 months from date of receipt.