

Product Datasheet

Prolactin peptide (orb372754)

Description	Prolactin peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Pro-Asp-Ile-Asn-Pro-Ala-Trp-Tyr-Ala-Ser-Arg-Gly-Ile-Arg-Pro-Val-Gly-Arg-Phe-NH ₂
Formula	C ₁₀₄ H ₁₅₈ N ₃₂ O ₂₆
Purity	> 95%
Source	Synthetic
MW	2272.62
CAS Number	235433-36-0
Expiration Date	12 months from date of receipt.

