

Product Datasheet

PRRS-RSAB-C (orb372773)

Description	Peptide to PRRS-RSAB-C
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Leu-Asn-Gly-Thr-Asp-Trp-Leu-Ser-Asn-His-Phe-Asp-Trp-Ala
Formula	$C_{77}H_{102}N_{20}O_{23}$
Purity	> 95%
Source	Synthetic
MW	1675.79
Expiration Date	12 months from date of receipt.