

Product Datasheet

pTH (3-34) (Bovine) (orb372778)

Description	Peptide to pTH (3-34) (Bovine)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Glu-Ile-Gln-Phe-Met-His-Asn-Leu-Gly-Lys-His-Leu-Ser-Ser-Met-Glu-Arg-Val-Glu-Trp-Leu-Arg-Lys-Lys-Leu-Gln-Asp-Val-His-Asn-Phe
Formula	$C_{175}H_{274}N_{52}O_{48}S_2$
Purity	> 95%
Source	Synthetic
MW	3938.55
Expiration Date	12 months from date of receipt.