

Product Datasheet

pTH-Related Protein Splice Isoform 3 (140-173) (orb372780)

Description	Peptide to pTH-Related Protein Splice Isoform 3 (140-173)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Ala-Leu-Leu-Trp-Gly-Leu-Lys-Lys- Lys-Lys-Glu-Asn-Asn-Arg-Arg-Thr- His-His-Met-Gln-Leu-Met-Ile-Ser-Leu- Phe-Lys-Ser-Pro-Leu-Leu-Leu-Leu
Formula	$C_{186}H_{313}N_{53}O_{44}S_2$
Purity	> 95%
Source	Synthetic
MW	4059.99
CAS Number	139872-85-8
Expiration Date	12 months from date of receipt.