

Product Datasheet

RES-701-1 (orb372785)

Description	Peptide to RES-701-1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Asn-Trp-His-Gly-Thr-Ala-Pro-Asp-Trp-Phe-Phe-Asn-Tyr-Tyr-Trp
Formula	$C_{103}H_{117}N_{23}O_{24}$
Purity	> 95%
Source	Synthetic
MW	2061.22
CAS Number	151308-34-8
Expiration Date	12 months from date of receipt.

