

## Product Datasheet

### RS domain derived peptide (orb372793)

|                         |                                              |
|-------------------------|----------------------------------------------|
| <b>Description</b>      | Peptide to RS domain derived peptide         |
| <b>Form/Appearance</b>  | Lyophilized powder                           |
| <b>Storage</b>          | Shipped at 4°C. Store at -20°C for one year. |
| <b>Note</b>             | For research use only                        |
| <b>Protein Sequence</b> | Gly-Arg-Ser-Arg-Ser-Arg-Ser-Arg-Ser-Arg      |
| <b>Formula</b>          | $C_{44}H_{85}N_{25}O_{15}$                   |
| <b>Purity</b>           | > 95%                                        |
| <b>Source</b>           | Synthetic                                    |
| <b>MW</b>               | 1204.33                                      |
| <b>Expiration Date</b>  | 12 months from date of receipt.              |