

Product Datasheet

Saposin C18 (orb372802)

Description	Peptide to Saposin C18
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Val-Lys-Glu-Val-Thr-Lys-Leu-Ile-Asp-Asn-Asn-Lys-Thr-Glu-Lys-Glu-Ile-Leu
Formula	$C_{93}H_{164}N_{24}O_{31}$
Purity	> 95%
Source	Synthetic
MW	2114.49
Expiration Date	12 months from date of receipt.

