

Product Datasheet

Saposin C22 (orb372803)

Description	Peptide to Saposin C22
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Cys-Glu-Phe-Leu-Val-Lys-Glu-Val-Thr-Lys-Leu-Ile-Asp-Asn-Asn-Lys-Thr-Glu-Lys-Glu-Ile-Leu
Formula	$C_{116}H_{196}N_{28}O_{37}S$
Purity	> 95%
Source	Synthetic
MW	2607.08
Expiration Date	12 months from date of receipt.

