

Product Datasheet

SH2 Domain Ligand (1) (orb372816)

Description	Peptide to SH2 Domain Ligand (1)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Asp-pTyr-Val-Pro-Met-Leu-NH ₂
Formula	C ₃₆ H ₅₆ N ₇ O ₁₃ PS
Purity	> 95%
Source	Synthetic
MW	857.93
Expiration Date	12 months from date of receipt.

