

Product Datasheet

SHNA (orb372820)

Description	Peptide to SHNA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Met-Ala-Pro-Arg-Gly-Phe-Ser-Ala-Leu-Leu-Leu-Leu-Thr-Ser-Glu-Ile-Asp-Leu-Pro-Val-Lys-Arg-Arg-Ala
Formula	$C_{119}H_{204}N_{34}O_{32}S$
Purity	> 95%
Source	Synthetic
MW	2655.23
Expiration Date	12 months from date of receipt.