

Product Datasheet

Sodefrin (orb372829)

Description	Peptide to Sodefrin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Ile-Pro-Ser-Lys-Asp-Ala-Leu-Leu-Lys
Formula	$C_{48}H_{86}N_{12}O_{15}$
Purity	> 95%
Source	Synthetic
MW	1071.28
CAS Number	161928-96-7
Expiration Date	12 months from date of receipt.