

Product Datasheet

Substance P-Gly-Lys-Arg (orb372853)

Description	Peptide to Substance P-Gly-Lys-Arg
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Pro-Lys-Pro-Gln-Gln-Phe-Phe-Gly-Leu-Met-Gly-Lys-Arg
Formula	$C_{77}H_{124}N_{24}O_{17}S$
Purity	> 95%
Source	Synthetic
MW	1690.06
Expiration Date	12 months from date of receipt.

