

Product Datasheet

SV40 Nuclear Transport Signal Peptide Analog (orb372858)

Description	Peptide to SV40 Nuclear Transport Signal Peptide Analog
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Cys-Gly-Tyr-Gly-Pro-Lys-Lys-Lys-Arg- Lys-Val-Gly-Gly
Formula	$C_{60}H_{104}N_{20}O_{15}S$
Purity	> 95%
Source	Synthetic
MW	1377.68
CAS Number	104914-40-1
Expiration Date	12 months from date of receipt.

