

Product Datasheet

Terlipressin Acetate (orb372865)

Description	Peptide to Terlipressin Acetate
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Gly-Gly-c[Cys-Tyr-Phe-Gln-Asn-Cys]-Pro-Lys-Gly-NH ₂
Formula	C ₅₂ H ₇₄ N ₁₆ O ₁₅ S ₂
Purity	> 95%
Source	Synthetic
MW	1227.37
CAS Number	14636-12-5
Expiration Date	12 months from date of receipt.