

Product Datasheet

Tetanus Toxin Peptide (orb372867)

Description	Peptide to Tetanus Toxin Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gln-Tyr-Ile-Lys-Ala-Asn-Ser-Lys-Phe-Ile-Gly-Ile-Phe-Glu
Formula	$C_{79}H_{120}N_{18}O_{21}$
Purity	> 95%
Source	Synthetic
MW	1657.95
Expiration Date	12 months from date of receipt.