

Product Datasheet

TNF-Alpha (10-36) (orb372871)

Description	Peptide to TNF-Alpha (10-36)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Lys-Pro-Val-Ala-His-Val-Val-Ala-Asn-Pro-Gln-Ala-Glu-Gly-Gln-Leu-Gln-Trp-Leu-Asn-Arg-Arg-Ala-Asn-Ala-Leu
Formula	$C_{131}H_{211}N_{43}O_{38}$
Purity	> 95%
Source	Synthetic
MW	2996.41
Expiration Date	12 months from date of receipt.

