

Product Datasheet

TNF-Alpha (71-82) (orb372874)

Description	Peptide to TNF-Alpha (71-82)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Pro-Leu-Ala-Gln-Ala-Val-Arg-Ser-Ser-Ser-Arg
Formula	$C_{51}H_{91}N_{19}O_{18}$
Purity	> 95%
Source	Synthetic
MW	1258.41
Expiration Date	12 months from date of receipt.

