

Product Datasheet

Tyr-Leptin (26-39) (orb372891)

Description	Peptide to Tyr-Leptin (26-39)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Lys-Val-Gln-Asp-Asp-Thr-Lys-Thr-Leu-Ile-Lys-Thr-Ile-Val
Formula	$C_{80}H_{137}N_{19}O_{25}$
Purity	> 95%
Source	Synthetic
MW	1765.10
Expiration Date	12 months from date of receipt.

