

Product Datasheet

Tyrosine Kinase Peptide 3, Acetylated, Amide, Phosphorylated (orb372895)

Description	Peptide to Tyrosine Kinase Peptide 3, Acetylated, Amide, Phosphorylated
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Arg-Arg-Leu-Ile-Glu-Asp-Ala-Glu-pTyr-Ala-Ala-Arg-Gly-NH ₂
Formula	C ₆₆ H ₁₁₀ N ₂₃ O ₂₄ P
Purity	> 95%
Source	Synthetic
MW	1640.75
Expiration Date	12 months from date of receipt.

