

Product Datasheet

Tyr-Proinsulin C-peptide (55-89) (orb372898)

Description	Peptide to Tyr-Proinsulin C-peptide (55-89)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Arg-Arg-Glu-Ala-Glu-Asp-Leu-Gln-Val-Gly-Gln-Val-Glu-Leu-Gly-Gly-Gly-Pro-Gly-Ala-Gly-Ser-Leu-Gln-Pro-Leu-Ala-Leu-Glu-Gly-Ser-Leu-Gln-Lys-Arg
Formula	$C_{162}H_{268}N_{50}O_{54}$
Purity	> 95%
Source	Synthetic
MW	3780.24
CAS Number	139532-11-9
Expiration Date	12 months from date of receipt.