

Product Datasheet

V5 Epitope Tag (orb372907)

Description	Peptide to V5 Epitope Tag
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Lys-Pro-Ile-Pro-Asn-Pro-Leu-Leu-Gly-Leu-Asp-Ser-Thr
Formula	$C_{64}H_{1081}N_{16}O_{20}$
Purity	> 95%
Source	Synthetic
MW	1421.67
Expiration Date	12 months from date of receipt.