

Product Datasheet

Vasopressin peptide (orb372909)

Description	Peptide to Val-Asp-(Arg8)-Vasopressin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Val-Asp-Cys-Tyr-Phe-Gln-Asn-Cys-Pro-Arg-Gly-NH ₂ (Disulfide bridge Cys3-Cys8)
Formula	C ₅₅ H ₇₉ N ₁₇ O ₁₆ S ₂
Purity	> 95%
Source	Synthetic
MW	1298.48
CAS Number	100930-18-5
Expiration Date	12 months from date of receipt.