

Product Datasheet

Valosin Peptide (orb372911)

Description	Peptide to Valosin Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Val-Gln-Tyr-Pro-Val-Glu-His-Pro-Asp-Lys-Phe-Leu-Lys-Phe-Gly-Met-Thr-Pro-Ser-Lys-Gly-Val-Leu-Phe-Tyr
Formula	$C_{141}H_{207}N_{31}O_{35}S$
Purity	> 95%
Source	Synthetic
MW	2928.50
CAS Number	99510-37-9
Expiration Date	12 months from date of receipt.

