

Product Datasheet

VIP Receptor-Binding Inhibitor L-8-K (orb372917)

Description	Peptide to VIP Receptor-Binding Inhibitor L-8-K
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Leu-Met-Tyr-Pro-Thr-Tyr-Leu-Lys
Formula	$C_{50}H_{77}N_9O_{12}S$
Purity	> 95%
Source	Synthetic
MW	1028.29
CAS Number	120550-85-8
Expiration Date	12 months from date of receipt.

