

Product Datasheet

VIP-Lys (Biotin) (orb372918)

Description	Peptide to VIP-Lys (Biotin)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Ser-Asp-Ala-Val-Phe-Thr-Asp-Asn-Tyr-Thr-Arg-Leu-Arg-Lys-Gln-Met-Ala-Val-Lys-Lys-Tyr-Leu-Asn-Ser-Ile-Leu-Asn-Lys (Biotin)
Formula	$C_{163}H_{263}N_{47}O_{46}S_2$
Purity	> 95%
Source	Synthetic
MW	3681.33
Expiration Date	12 months from date of receipt.