

Product Datasheet

Tyr-CRF (orb372935)

Description	Peptide to Tyr-CRF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Ser-Glu-Glu-Pro-Pro-Ile-Ser-Leu-Asp-Leu-Thr-Phe-His-Leu-Leu-Arg-Glu-Val-Leu-Glu-Met-Ala-Arg-Ala-Glu-Gln-Leu-Ala-Gln-Gln-Ala-His-Ser-Asn-Arg-Lys-Leu-Met-Glu-Ile-Ile-NH ₂
Formula	C ₂₁₇ H ₃₅₃ N ₆₁ O ₆₅ S ₂
Purity	> 95%
Source	Synthetic
MW	4920.73
CAS Number	100513-58-4
Expiration Date	12 months from date of receipt.