

Product Datasheet

Beta-Amyloid (1-43) (orb372937)

Description	Peptide to Beta-Amyloid (1-43)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Ala-Glu-Phe-Arg-His-Asp-Ser-Gly-Tyr-Glu-Val-His-His-Gln-Lys-Leu-Val-Phe-Phe-Ala-Glu-Asp-Val-Gly-Ser-Asn-Lys-Gly-Ala-Ile-Ile-Gly-Leu-Met-Val-Gly-Gly-Val-Val-Ile-Ala-Thr
Formula	$C_{207}H_{318}N_{56}O_{62}S$
Purity	> 95%
Source	Synthetic
MW	4615.24
CAS Number	134500-80-4
Expiration Date	12 months from date of receipt.