

Product Datasheet

Parathyroid Hormone (1-44) (orb372944)

Description	Peptide to Parathyroid Hormone (1-44)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Val-Ser-Glu-Ile-Gln-Leu-Met-His-Asn-Leu-Gly-Lys-His-Leu-Asn-Ser-Met-Glu-Arg-Val-Glu-Trp-Leu-Arg-Lys-Lys-Leu-Gln-Asp-Val-His-Asn-Phe-Val-Ala-Leu-Gly-Ala-Pro-Leu-Ala-Pro-Arg
Formula	$C_{225}H_{366}N_{68}O_{61}S_2$
Purity	> 95%
Source	Synthetic
MW	5063.93
CAS Number	85568-24-7
Expiration Date	12 months from date of receipt.