

Product Datasheet

FDC-SP (30-85) (orb372956)

Description	Peptide to FDC-SP (30-85)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Ile-Ser-Asp-Ser-Asp-Glu-Leu-Ala-Ser-Gly-Phe-Phe-Val-Phe-Pro-Tyr-Pro-Tyr-Pro-Phe-Arg-Pro-Leu-Pro-Pro-Ile-Pro-Phe-Pro-Arg-Phe-Pro-Trp-Phe-Arg-Arg-Asn-Phe-Pro-Ile-Pro-Ile-Pro-Glu-Ser-Ala-Pro-Thr-Thr-Pro-Leu-Pro-Ser-Glu-Lys
Formula	$C_{314}H_{447}N_{71}O_{78}$
Purity	> 95%
Source	Synthetic
MW	6464.49
Expiration Date	12 months from date of receipt.