

Product Datasheet

I-309 (orb372964)

Description	Peptide to I-309
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Ser-Met-Gln-Val-Pro-Phe-Ser-Arg-Cys-Cys-Phe-Ser-Phe-Ala-Glu-Gln-Glu-Ile-Pro-Leu-Arg-Ala-Ile-Leu-Cys-Tyr-Arg-Asn-Thr-Ser-Ser-Ile-Cys-Ser-Asn-Glu-Gly-Leu-Ile-Phe-Lys-Leu-Lys-Arg-Gly-Lys-Glu-Ala-Cys-Ala-Leu-Asp-Thr-Val-Gly-Trp-Val-Gln-Arg-His-Arg-Lys-Met-Leu-Arg-His-Cys-Pro-Ser-Lys-Arg-Lys
Formula	$C_{370}H_{604}N_{114}O_{99}S_8$
Purity	> 95%
Source	Synthetic
MW	8490.11
Expiration Date	12 months from date of receipt.