

Product Datasheet

Tyr-Uroguanylin (orb372985)

Description	Peptide to Tyr-Uroguanylin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Thr-Asp-Glu-Cys-Glu-Leu-Cys-Ile-Asn-Val-Ala-Cys-Thr-Gly-Cys (Disulfide bridge Cys4-Cys12, Cys7-Cys15)
Formula	$C_{69}H_{105}N_{17}O_{27}S_4$
Purity	> 95%
Source	Synthetic
MW	1732.96
Expiration Date	12 months from date of receipt.

