


Product Datasheet

SLC12A6 antibody (Biotin) (orb372992)

Description	Goat polyclonal antibody against SLC12A6 conjugated to Biotin
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	KCC3 / SLC12A6, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 120kDa band observed in Human Heart lysates (calculated MW ranging 121-128kDa according to NP_598408.1; NP_005126.1; NP_001035959.1; NP_001035960.1; NP_001035961.1; NP_001035962.1). See non-biotinylated parental product's datasheet for further Q
Protein Sequence	SRHIDVCSKTKE.
Clonality	Polyclonal
MW	128; 122; 121; 127; 126
Entrez	9990
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...