


Product Datasheet

SNAP25 antibody (Biotin) (orb372995)

Description	Goat polyclonal antibody against SNAP25 conjugated to Biotin
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SNAP25, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 26kDa band observed in Rat Brain lysates (calculated MW of 23.3kDa according to NP_003072.2 and NP_570824.1). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	DEANQRATKMLGSG.
Clonality	Polyclonal
MW	23.3; 23.3; 23.3
Entrez	6616
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Brain lysat...