


Product Datasheet

PDIA2 antibody (Biotin) (orb373000)

Description	Goat polyclonal antibody against PDIA2 conjugated to Biotin
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PDIA2 (aa477-481), Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 65kDa band observed in Rat Pancreas lysates (calculated MW of 58.5kDa according to Rat NP_001099245.2). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.03-0.1µg/ml.
Protein Sequence	EYKSTRDLETFSK.
Clonality	Polyclonal
MW	58.2
Entrez	64714
Dilution Range	ELISA: 1:64000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Pancreas ly...