



Product Datasheet

PARP1 antibody (orb373014)


Description	Rabbit polyclonal antibody to PARP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human PARP1(NP_001609.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	89KDa/113KDa
Uniprot ID	P09874
Dilution Range	WB: 1:1000-3000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various cell lines.



Immunohistochemistry of paraffin-embedded tissue sections.



Immunohistochemistry of paraffin-embedded tissue sections.