


Product Datasheet

RTN4RL1 antibody (orb373021)

Description	Rabbit polyclonal antibody to RTN4RL1
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human RTN4RL1(NP_848663.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49kDa
Uniprot ID	Q86UN2
Dilution Range	WB: 1:1000-4000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...