


Product Datasheet

LILRB2 antibody (orb373063)

Description	Rabbit polyclonal antibody to LILRB2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human LILRB2/CD85d/ILT4(Q8N423).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	Q8N423
Dilution Range	WB: 1:1000-4000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...