


## Product Datasheet

### RBM34 antibody (orb373065)

|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to RBM34                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                       |
| <b>Conjugation</b>         | Unconjugated                                            |
| <b>Tested Applications</b> | IHC, WB                                                 |
| <b>Immunogen</b>           | Recombinant fusion protein of human RBM34(NP_055829.2). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.      |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.             |
| <b>Note</b>                | For research use only                                   |
| <b>Isotype</b>             | IgG                                                     |
| <b>Clonality</b>           | Polyclonal                                              |
| <b>MW</b>                  | 48kDa                                                   |
| <b>Uniprot ID</b>          | <b>P42696</b>                                           |
| <b>Dilution Range</b>      | WB: 1:500-2000                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                         |



Western blot analysis of extracts of var...