



Product Datasheet

DYNC111 antibody (orb373160)


Description	Rabbit polyclonal antibody to DYNC111
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human DYNC111(NP_001129028.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	73kDa
Uniprot ID	O14576
Dilution Range	WB: 1:200-2000, IHC-P: 1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extract of vari...



Immunohistochemical staining of human ga...



Immunohistochemical staining of rat brai...