




Product Datasheet

AKT2 antibody (orb37319)

Description	Rabbit polyclonal antibody to AKT2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	AKT2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55769
Uniprot ID	P31751
NCBI	NP_001617.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, FACS: 1:10-50, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of AKT2 antibody



Immunofluorescence analysis of 293 cell

...