


Product Datasheet

EIF4A3 antibody (orb373215)

Description	Rabbit polyclonal antibody to EIF4A3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human eIF4A3(NP_055555.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47kDa
Uniprot ID	P38919
Dilution Range	WB: 1:1000-3000, IHC-P: 1:50-200, IF/ICC: 1:50-100, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...