



Product Datasheet

USH1C antibody (orb373218)


Description	Rabbit polyclonal antibody to USH1C
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of zebrafish USH1C(NP_001035018.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	73kDa
Dilution Range	WB: 1:200-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat hear...



Immunohistochemical staining of mouse ki...



Immunohistochemical staining of human st...