



Product Datasheet

SAP18 antibody (orb373222)


Description	Rabbit polyclonal antibody to SAP18
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human SAP18(NP_005861.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18kDa
Uniprot ID	O00422
Dilution Range	WB: 1:200-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extract of vari...



Immunohistochemical staining of rat brai...



Immunohistochemical staining of mouse br...