


Product Datasheet

NCAPH antibody (orb373244)

Description	Rabbit polyclonal antibody to NCAPH
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human NCAPH(NP_056156.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85kDa
Uniprot ID	Q15003
Dilution Range	WB: 1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extract of vari...