



Product Datasheet

COTL1 antibody (orb373245)


Description	Rabbit polyclonal antibody to COTL1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human COTL1(NP_066972.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	Q14019
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extract of vari...



Immunohistochemical staining of mouse br...



Immunohistochemical staining of rat kidn...